ABBREVIATION

- > Mm ~ mili meter
- ≻ Cm ~ Centimeter
- ➤ G Gram
- ≻ Mg mili gram
- ➢ µl micro liter
- ➤ mM ~ mili molar
- > F Bacterial strain isolated from fermented fish
- > KH Bacterial strain isolated from fermented mustard (kharoli)
- > RM Bacterial strain isolated from rice beer starter culture
- > D Bacterial strain isolated from curd (dahi)
- > M Bacterial strain isolated from milk

CONTENTS

CHAPTER-1

- 1.1 INTRODUCTION
- 1.2 Methods of classification
- 1.3 Databases

CHAPTER-2 Review & literature CHAPTER-3 Objective

CHAPTER-4 Materials and method

- 4.1 MATERIALS
- 4.1.1 Tested indicator strains
- 4.1.2 Fermented food
- 4.2 METHODS
- 4.2.1 Isolation of bacteria from sample by Serial dilution
- **4.2.2** Screening of isolated bacteria which shows antimicrobial activity
- **4.2.3** Determination of zone of inhibition through well diffusion assay:

4.2.4 Supernatant treated with protease or catalase to find out whether antimicrobial agent is protein or peroxide

- 4.2.5 EFFECT OF pH ON ANTIMICROBIAL ACTIVITY:
- 4.2.6 Effect of concentration on antimicrobial activity:
- 4.2.7 Characterization of bacteria by

CHAPTER 5 RESULT

5.1.1 Isolation of bacteria from curd and milk samples & antimicrobial assay

5.1.2 Seven & six bacterial strains were Isolated respectively from fermented bamboo shoot & soybean respectively and checked their antimicrobial activity.

5.1.3 Nine bacterial strains isolated from rice beer starter culture & check antimicrobial activity

5.1.4 Eleven bacterial strains were isolated from fermented fish samples & measure zone of inhibition

5 1.5 23 Bacterial strains isolated from fermented mustard & measured zone of inhibition.

5.2 An antimicrobial assay after treating with 20 μ L catalase (50U/ml) and trypsin (0.6mg/ml)

- 5.3 An effect of concentration on antimicrobial activity
- 5.4 An effect of ph on antimicrobial activity
- 5.5 Gram staining of different bacterial strain
- 5.6 Biochemical test
- 5.7 Biochemical assay of different bacterial strain:
- 5.8 Estimate size of protein by SDS PAGE.